ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE:

The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS:

Provides support to prevent further injury by stabilizing the ankle. It may be suitable for use following post-operative or acute injuries such as ankle sprains and Chronic Ankle Instability (CAI). Soft-good/semi-rigid designed to restrict motion through elastic or semi-rigid construction.

CONTRAINDICATIONS: N/A

WARNINGS AND PRECAUTIONS:

- Use only on ankle.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- Do not use if shells or straps are damaged.
- Do not use over open wounds.
- Do not use this device if it was damaged and/or package has been opened.
- If you develop an allergic reaction and/or experience itchy, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS:

Put on an absorbent sock

- Adjust Heelpad Width:**
 - Open brace straps. Lift bottom of aircells away from the shell to expose the heelpad.
 - Adjust heelpad width for a snug fit (Fig. 1) at the ankle. Firmly press heel straps and aircells back in place.
- Apply Brace:**
 - Step into the brace inline with the foot diagram on the heelpad (a sock or elastic wrap can be worn for additional comfort and compression).
 - Align brace sides with ankle (Fig. 2).
- Secure Straps:**
 - Secure straps from bottom to top (Fig. 3). Do not over tighten straps.
 - Put on laced shoe or sneaker.
- Adjust Brace:**
 - Squeeze brace sides together with one hand.
 - Re-Tighten straps from bottom to top with other hand (Fig. 3) until comfortable.

ANKLE BRACE FIT ADJUSTMENT: Ankle brace aircells are preinflated and normally do not require adjustment. If you experience pinching or uneven pressure, re-adjust straps and heelpad for proper fit. If condition persists, adjust aircell volume (**too much air will reduce support**):

- Unfold valve on top of aircell. (Fig. 4)
- Insert inflation tube 15 inches into valve.
- Remove a small amount of air by gently squeezing aircell or. Add a small amount of air by blowing into tube.
- Squeeze valve flat just below tube, and slowly remove tube. Pinch valve shut when removing the tube.
- Fold valve over and tuck against the shell.

HIGH ALTITUDE: At high altitudes the aircells will expand beyond their optimal level and air may have to be removed from the aircells (see ankle brace fit adjustment). When flying, re-adjust straps to a comfortable pressure.

CLEANING INSTRUCTIONS:

The entire ankle brace can be hand washed in lukewarm water with mild soap and air dried.

FOR SINGLE PATIENT USE ONLY.

~~NOT MADE WITH NATURAL RUBBER LATEX.~~

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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